



# COONEY ORTHODONTICS

## **Using the device:**

When used for 20 minutes per day, the AcceleDent Aura device will help your teeth respond to orthodontic treatment more quickly. It emits a gentle vibration with a stutter pulse at 10 minutes to let you know you're 1/2 way through the daily session. Try to find a time that you can use the device uninterrupted – it is best not to stop use in the middle of a session.

While wearing your Invisalign aligners, place the appliance in your mouth and bite down lightly, with just enough force to hold the appliance in place. The mouthpiece can become damaged by clenching it while using the AcceleDent. The appliance does not need to be held by hand once it's placed inside the mouth. You will be more comfortable if you relax the jaw and facial muscles, while maintaining enough force to stabilize the mouthpiece. At the end of the session, detach the mouthpiece, making sure NOT to rock it side to side or up and down. Rinse the mouthpiece in cool water and gently wipe it with a soft cloth if needed. It is important to keep the activator handle dry; only wipe it with a soft, dry cloth.


## **Charging and Powering the Device:**

Fully charge the Activator prior to initial use by charging it overnight. The Activator is properly charging when it is plugged into the Power Adapter and the LED Indicator is slowly flashing orange. To maximize the time between charges, the Activator should be recharged every 3 to 5 days or when the LED Indicator is quickly flashing orange when the Activator is powered on but not charging. Do not let the battery completely discharge as this may affect the data reporting of the FastTrac Usage Report.

***DO NOT USE THE ACCELEDEMENT AURA WHEN THE ACTIVATOR IS CHARGING OR WHEN IT IS CONNECTED TO THE COMPUTER***

## **Contraindications for Use**

- Use of osteoporosis drugs
- Poor oral hygiene
- Periodontal disease that is not under full control for at least 3-4 months prior to the start of treatment



# COONEY ORTHODONTICS

## **AcceleDent Aura Agreement:**

Dr. Cooney has prescribed the use of the AcceleDent Aura appliance to help your Orthodontic treatment work faster and provide a better experience from start to finish. It is important to maintain a regular schedule of 20 minute daily use and proper care and storage of the appliance.

If for some reason, the AcceleDent Aura has been broken and is not covered under the manufacturer's warranty or the device is lost, a \$700 replacement fee will apply.

If for some reason, the mouthpiece is lost or damaged, a \$150 replacement fee will apply.

All warranty information is enclosed in packaging.

I have read and understood the AcceleDent Aura agreement.

---

Date

---

Patient Name

---

Patient Signature

---

Parent Signature